

# PERFIT

**DA DENTAL  
ADVISOR**  
Product insights you can trust.

## Hydrophilic Impression Materials



With over 10 years of clinical success, PERFIT has become one of the most reliable and widely used impression materials on the market. HUGE continues to expand their range to meet the diverse needs of dental professionals — ensuring there's always a [PERFIT solution for every indication](#).

### Outstanding Hydrophilicity & Thixotropy

PERFIT offers excellent flow in moist environments, enabling precise detail reproduction. Its thixotropic properties prevent dripping, making the procedure more comfortable for both dentists and patients.

### Excellent dimensional stability

Low shrinkage and high dimensional stability allow impressions to remain accurate for up to **28 days** - ideal for both chairside efficiency and lab processing.

Elastic recovery > 99.5%

Linear dimensional change < 0.2%

### Smart Handling & Reliable Results

Thanks to its elastic memory, PERFIT can be easily removed from the mouth and returns to its original shape. Available in multiple viscosities and hardness levels, PERFIT provides the flexibility clinicians need to achieve high-precision impressions.

### Doctor/Patient-Friendly Working & Setting Times

PERFIT offers efficient curing with sufficient working time and a shorter intraoral setting time. This makes procedures easier for dentists and more comfortable for patients. Faster cross-linking also helps reduce inaccuracies caused by micro-movement.

### Optimized Timing for Every Impression Technique

PERFIT Viscosity	Putty		Heavy Body			Regular Body (Transparent)	Light Body		Super Light Body	Bite Registration	
	Type 0		Type 1			Type 2	Type 3		Type 3	Bite	
Working Time	Normal Set	Fast Set	Dispenser-mix Normal Set	Dispenser-mix Fast Set	Auto-mix Fast Set	Fast Set	Normal Set	Fast Set	Normal Set	Normal Set	Fast Set
	1:30 min.	1:15 min.	2:00 min.	1:30 min.	1:30 min.	1:30 min.	2:15 min.	1:30 min.	2:15 min.	0:30 min.	0:12 min.
Setting Time in Mouth	3:00 min.	2:30 min.	3:30 min.	3:00 min.	3:00 min.	3:00 min.	3:30 min.	2:30 min.	3:30 min.	1:30 min.	0:45 min.
Package	400 g base + 400 g catalyst + 2 measuring spoons		50 mL × 2 cartridges 50 mL × 2 cartridges + 6 mixers			380 mL × 1 cartridge	50 mL × 2 cartridges	50 mL × 2 cartridges 50 mL × 2 cartridges + 12 mixers + 12 oral tips		50 mL × 2 cartridges 50 mL × 2 cartridges + 12 mixers + 12 oral tips	

\*Data from HUGE Laboratory

# HUGE

HUGE DENTAL USA LLC

Web: [www.hugedentalusa.com](http://www.hugedentalusa.com)

Add: 4010 Valley Blvd STE 104, Walnut, CA 91789, USA

Email: [info@hugedentalusa.com](mailto:info@hugedentalusa.com) Tel: +1 626 283 5808



Follow @HUGE Dental USA

v2.0-202507