

PERFIT

Elastomeric
Impression Material

MONOPHASE

Type 2 Regular Body Impression Material



Precision and Reliability Combined in Every Single Impression

Would you like to simplify impression making while ensuring outstanding accuracy?

PERFIT Monophase is the perfect choice. With its single-viscosity formulation, it enables one-step impressions, saving time without compromising precision.

Its excellent hydrophilicity, high tear resistance, and dimensional stability provide reliable results for *implant, crown, and bridge impressions*— every time.

PERFIT MONOPHASE

INDICATIONS

- ✓ For crown, bridge, inlay, onlay and implant impressions
- ✓ Recommended for one-step impression technique

BENEFITS


- ✓ High hydrophilicity and thixotropy
- ✓ Moderate flow for high impression accuracy
- ✓ Excellent handling properties, easy to use
- ✓ High tear and tensile strength
- ✓ Good dimensional stability - up to 28 days

TECHNICAL PARAMETERS

Product	Monophase-Normal Set
Working Time	2min00s
Setting Time in Oral Cavity	3min30s
Elastic Recovery	> 99.5%
Dimensional Change Rate	< 0.2%
Hardness Shore A(24 h)	≈60
Color	Light Yellow
Shelf Life	3 Years

*Data from VinciSmile Lab.

Key Mechanical Properties Comparison



OPERATION GUIDE



Working Time
2'00"

Minimum Time In
The Oral Cavity 3'30"

PACKAGING

Standard: 50 mL × 2
cartridges
Kit: 50 mL × 2 cartridges
+ 6 mixers

Single Implant Impression

Multiple / Full-Arch Implant Impression

The result: Tear-resistant and elastic.

PERFIT Monophase offers excellent elastic recovery and tear resistance, making impressions easy to remove from the mouth without distortion. Its optimal hardness also provides stability for transfer posts in implant impressions, ensuring accurate details and reliable final results.

HUGE

HUGE DENTAL USA LLC

Web: www.hugedentalusa.com

Add: 4010 Valley Blvd STE 104, Walnut, CA 91789, USA

Email: info@hugedentalusa.com Tel: +1 626 283 5808



Follow @HUGE Dental USA

v2.0-202507