

Dental Model

Instructions for Use

Version: 04
Issuance Date: Nov. 20, 2025

This product should be used by the operator in accordance with the instructions of the user manual and within the scope of application. The operator should read this manual carefully before using it.

【Product Name】 Dental Model

【Structure and composition】 This product is a resin material.

【Intended purpose】 Used to make dental models

【Model】 M 2.0, M 3.0, H 2.0, D 2.0, L 2.0

【Operation description】

1. This product is a 3D printing material suitable for UV light curable 3D printers. Please shake well before use.
2. Please wear goggles and rubber gloves when using.
3. The cleaning of L 2.0 models is performed with fresh water instead of alcohol.

【Storage conditions and methods】

The product should be kept away from light, moisture and sealed under the condition of 5-25 C.

【Precautions】

Those who are allergic to acrylic resin ingredients should not be allowed. If allergic reactions occur during use, please consult a doctor.


Notification of any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority in your country.

【Adverse events】 There are no clinical cases of adverse events for this product.

【Shelf life】 12 months

【Accessories】 None

【Symbols used in packaging/labeling】



 Shandong Huge Dental Material Corporation
No. 68 Shanhai Road, Donggang District, Rizhao City, Shandong Province, 276800, P.R. China
Tel: +86 633 2277268
www.hugedental.com marketing@hugedental.com

 MedNet EC-REP GmbH
Borkstrasse 10, 48163 Muenster, Germany

